

## MEP Instruments Media Release

### For immediate publication

Media Contact:

Daniel Ledergerber, Creative Spiders Sydney, 02 9590 4777, [ledi@creativespiders.com](mailto:ledi@creativespiders.com)

Technical Contact Metrohm Titration:

Carlo Chiara, MEP Instruments Victoria, 03 9810 4516, [cc@mep-instruments.com.au](mailto:cc@mep-instruments.com.au)

Photo: Metrohm Titroprocessor 796

Number of words: 312

Media release can be downloaded from <http://www.mep.net.au/press/release>

### Metrohm Potentiometric Titration



#### **796 Titroprocessor with Metrohm Sample changer for state-of-the art laboratory automation**

The 796 Titroprocessor is Metrohm's easy-to-use top-class titration and measuring instrument for analyses in the laboratory and plant. This instrument makes routine operation extremely simple. Thanks to its efficiently used windows, hardkeys and softkeys routine tasks can be carried out after a minimal familiarization period and without the time-consuming study of operating instructions. The instrument can of course also be run in expert mode with access to all parameters. A personal key card provides access to the 796 Titroprocessor and all important functions. The instrument reads your preferred configuration from your key card and presents itself to you exactly as you want it. Methods created on the 726 Titroprocessor can be transferred directly to the 796 Titroprocessor.

Selected features:

- Key Card with your personal instrument settings
- Backlit LCD display - shows operating procedures, determination sequences, real-time curves and results
- Tactile touch keyboard
- 1 or 2 (option) galvanically separated measurement groups with various Titration and measurement modes: DET, MET, SET, SEC, KFT, KFC, MEAS, CAL
- Designed for GLP: Service interval monitoring and GLP test for pH electrodes
- Optional built-in printer or commercial printer

MEP Instruments is sole distributor of leading names in sample analysis technology in Australia and New Zealand. The joint venture company formed by Metrohm Switzerland



MEP Instruments Media Release p. 2

and Anton Paar Austria provides comprehensive service, support and training for the whole product line of Metrohm, Anton Paar, Autolab, Camag, Ismatec, Physica and Retsch in Australia and New Zealand. Contact MEP in Australia: MEP Instruments Pty Ltd (Australia) Unit 5, 112-118 Talavera Rd, North Ryde, NSW 2113 Tel 02 9878 6900 Fax 02 9878 6999 [info@mep-instruments.com.au](mailto:info@mep-instruments.com.au); or in New Zealand: New Zealand Brooke House, 24 Balfour Road, Parnell, Auckland Tel (+64) 09 366 1236 Fax (+64) 09 366 1235 Email [info@mep-instruments.co.nz](mailto:info@mep-instruments.co.nz).  
Visit Mep online: <http://www.mep.net.au>