

BLACKOUT PERIOD ENDS 20/09/2011

## Metrohm Titration

Introducing the Ti-Touch: Maximum simplicity and minimal footprint

Metrohm Ltd, the Swiss company renowned for its reliable ion analysis systems, is introducing its latest compact titrators. "Titration Reduced to the Max" is the chosen inauguration slogan. The resourceful engineers at Metrohm have indeed managed to reduce size, clutter and handling complexity to an absolute minimum – thus maximal reduction!

Buret, stirrer, dosing and control modules are all contained in a single unit. The Ti-Touch connects directly with your network, no need for a PC. It saves forge-proof PDF reports directly on a USB flash drive or in the intranet.

Patented contact-free reagent exchange guarantees user safety. A single touch on one of 14 icons on the screen starts the selected process ensuring maximum simplicity.

The 916 Ti-Touch is the compact titrator for potentiometric titration incorporating 4 different titration modes. The new compact titrator comes also in the versions 916 Salt Ti-Touch, 916 Oil Ti-Touch and 916 Food Ti-Touch, which are all-inclusive packages containing the titrator plus the complete accessories for the respective applications.

The 915 KF Ti-Touch is the version for Karl Fischer Titration. As the Karl Fischer titration uses toxic solutions, the 915 KF Ti-Touch lays the greatest emphasis on user safety. Special KF icons inform whether the instrument is still being conditioned or whether the water determination can already be started. The 915 can administer and communicate data just as well as the variant for potentiometric titration. With the Ethernet interface, the 915 KF Ti-Touch incorporates easily into a network or connects to a LIMS. Peripheral instruments connect via the built-in USB port. The Ti-Touch is available in Australia and New Zealand at MEP Instruments.

### **MEP Instruments Pty Ltd**

PO Box 2173, Boronia Park  
Unit 11, 56 Buffalo Road  
Gladesville, NSW 2111  
**T** (+61) 02 8899 5200  
**F** (+61) 02 8899 5299  
**E** info@mep.net.au  
**www.mep.net.au**

**ABN** 93 081 861 645

**MEP Instruments Pty Ltd**

PO Box 2173, Boronia Park  
Unit 11, 56 Buffalo Road  
Gladesville, NSW 2111  
**T** (+61) 02 8899 5200  
**F** (+61) 02 8899 5299  
**E** [info@mep.net.au](mailto:info@mep.net.au)  
**www.mep.net.au**

**ABN** 93 081 861 645

---

Available in Australia and New Zealand at MEP Instruments.

For further information please contact MEP Instruments in Australia or in New Zealand:

AU:

MEP Instruments Pty Ltd , Unit 11+12 , 56 Buffalo Road , Gladesville NSW  
2111 AUSTRALIA , Phone: +61 2 9878 6900 , Fax : +61 2 9878 6999,  
[info@mep.net.au](mailto:info@mep.net.au), [www.mep.net.au](http://www.mep.net.au)

NZ:

MEP Instruments Ltd (New Zealand) Building 2, 5 Ceres Court Mairangi Bay North  
Shore City 0630 Tel 09 912 1330 Fax 09 912 6892 Email [info@mep.net.nz](mailto:info@mep.net.nz)  
[www.mep.net.nz](http://www.mep.net.nz)

Image 3376

Compact in size the new titrators from Metrohm fit in any lab and are nice to their operators.

Image 3365

Maximal connectivity and integration into networks with multiple ports.