

MEP Instruments Media Release

“pH Measurement – Buffer Solutions in Sachets”

For immediate publication

Number of words: 368

Media & Technical Contact:

Carlo Chiara cc@mep.net.au

Text and pictures are available on our website: <http://mep.net.au/pressreleases.php>

Photos:



Metrohm Buffer solutions pH 4.00 / 7.00 / 9.00 in disposable sachets

Company Logos: MEP / Metrohm



pH measurement made even easier:

Metrohm Buffer solutions pH 4.00 / 7.00 / 9.00 in disposable sachets

The Metrohm «ready to use» buffer solutions pH 4.00, pH 7.00 and pH 9.00 are now available in disposable single sachets. The use of single packed buffer solutions guarantees that always fresh, contamination-free buffer solutions are used. The buffer solutions are passed through a sterile filter so that a shelf life of two years can be guaranteed.

The chemical composition is identical to the well established Metrohm buffer solutions 6.2307.1X0 in the 500 mL bottles. Naturally, the buffer solutions in the disposable sachets are fully traceable to NIST standards.

The following units are available:

- 6.2307.200 **Buffer solution pH 4.00** (25 °C), 30 x 30 mL,
shelf life 2 years, traceable to NIST
- 6.2307.210 **Buffer solution pH 7.00** (25 °C), 30 x 30 mL,
shelf life 2 years, traceable to NIST
- 6.2307.220 **Buffer solution pH 9.00** (25 °C), 30 x 30 mL,
shelf life 2 years, traceable to NIST
- 6.2307.230 **Buffer solutions pH 4.00 / 7.00 / 9.00**, 3 x 10 x 30
mL,
shelf life 2 years, traceable to NIST

For more information contact:

MEP Instruments Pty Ltd, PO Box 1880, North Ryde 2113, Phone 02 98786900. Fax: 02 98786999

MEP Instruments is sole distributor of leading names in sample analysis technology in Australia and New Zealand.

The joint venture company formed by Metrohm Switzerland and Anton Paar Austria provides COMPREHENSIVE service, support and training for the whole product line of Metrohm, Anton Paar, Applikon, Autolab, Camag, Ismatec, and Retsch in Australia and New Zealand.

Australia

MEP Instruments Pty Ltd (Australia), Unit 5, 112-118 Talavera Road, North Ryde, NSW 2113 Tel 02 9878 6900 Fax 02 9878 6999 info@mep.net.au;

New Zealand:

MEP Instruments Ltd (New Zealand) Brooke House, 24 Balfour Road, Parnell, Auckland Tel (+64) 09 366 1236, Fax (+64) 09 366 1235 info@mep.net.nz

Visit MEP online: <http://www.mep.net.au>