

MEP Instruments Media Release

“ISMATEC - NEW LABORATORY PUMPING CATALOG AVAILABLE NOW”

For immediate publication

Expires:

Number of words: 193

Media & Technical Contact:

Reto Broger rb@mep.net.au – 02 98786900



Photos:

10% introductory offer on low-cost Ecoflow

Ismatec introduces the all new low-cost ECOFLOW for flow rates from **100 ml/min to 3800 ml/min.**

Ismatec customer benefit from an additional discount of 10% on all ECOFLOW orders placed by 31 Dec 2005.

Ismatec liquid transfer program 2006

Ismatec pumps are available for virtually all liquid transfer applications with controlled flow rates of 0.001 ml/min to 45'000 ml/min. The comprehensive catalogue is available at MEP Instruments.

Company Logos



ISMATEC NEW LABORATORY PUMPING CATALOG IS AVAILABLE NOW

The high accuracy and long service-life of Ismatec pumps are well known worldwide.

Ismatec's latest 84 page product catalogue is now available from MEP Instruments. It features a wealth of information for users wanting to transfer liquid at rates from as low as 0.001 ml/min to 45'000 ml/min. The catalogue is split into 3 sections – tubing, gear and rotary piston pumps. They are all associated to very unique applications and will meet the most challenging demands. The amount of channels per pump varies between 1 and 24.

Tubing pumps are immune to dry running, require no flushing and no cleaning. Due to Ismatec's unique convex rollers liquid is treated gently, an ideal feature for the pumping of sensitive cell media.

Gear pumps are the pump of choice for pulse free pumping and to overcome differential pressures of up to 8.7 bar.

Rotary piston pumps show extreme resistance to corrosive and aggressive chemicals. They are also extremely accurate and highly reproducible results are achieved.

Special offer: MEP is offering an introductory 10% discount on all Ecoflow orders until the end of the year 2005.

For your free catalogue contact:

MEP Instruments Pty Ltd, PO Box 1880, North Ryde 2113, Phone 02 98786900. Fax: 02 98786999

MEP Instruments is sole distributor of leading names in sample analysis technology in Australia and New Zealand.

The joint venture company formed by Metrohm Switzerland and Anton Paar Austria provides COMPREHENSIVE service, support and training for the whole product line of Metrohm, Anton Paar, Applikon, Autolab, Camag, Ismatec, and Retsch in Australia and New Zealand.

AUSTRALIA

MEP Instruments Pty Ltd (Australia), Unit 5, 112-118 Talavera Road, North Ryde, NSW 2113 Tel 02 9878 6900 Fax 02 9878 6999 info@mep.net.au;

New Zealand:

MEP Instruments Ltd (New Zealand) Brooke House, 24 Balfour Road, Parnell, Auckland Tel (+64) 09 366 1236, Fax (+64) 09 366 1235 info@mep.net.nz

Visit MEP online: <http://www.mep.net.au>