

## MEP Instruments Media Release

For immediate publication

Media Contact:

Sinikka D teindre, MEP Instruments Pty Ltd Sydney 02 9878 6900, [sd@mep-instruments.com.au](mailto:sd@mep-instruments.com.au)  
Daniel Ledergerber, Creatvie Spiders Sydney, 02 9590 4777, [ledi@creativespiders.com.au](mailto:ledi@creativespiders.com.au)

Technical Contact Metrohm Titration:

Carlo Chiara, MEP Instruments Pty Ltd Victoria, 03 9810 4516, [cc@mep-instruments.com.au](mailto:cc@mep-instruments.com.au)

Photo: Metrohm KFD Titrino 758

Number of words: 230

Media release (incl. pictures in print resolution) can be downloaded from  
<http://www.mep.net.au/press/release>

---

### 758 KFD Titrino — 'live' Karl Fischer

stands for more information in KF titration. The graphical display of the reagent volume as a function of time allows to draw conclusions concerning the reaction kinetics and to detect side reactions. The latest member of our Titrino family with programmable stirrer control combines the advantages of the classical volumetric Karl Fischer titration with the advances of digital technology.

The KFD Titrino can control two additional dispensing instruments for titrations (685 Dosimats or 700 Dosinos) as well as two 776 or 765 auxiliary Dosimats. With one single KF titrator it is thus possible to have Karl Fischer reagents of different concentration or type (one-component or two component reagents) permanently available.



The 758 KFD Titrino affords the possibility to save your titration methods not only in the instrument's memory, but also on a memory card, thus making methods transportable.

The KFD Titrino offers extremely simple operation with the highest degree of comfort in the language of your choice: English, German, French, Spanish, Italian, Portuguese or Swedish.

Two RS 232C interfaces are available for the connection of a balance, a printer or a PC.

The complete automation of the Karl Fischer titration can be achieved using the 760 Titration Sample Changer and the Metrodata TiNet programme package. (Options).

Metrodata VESUV 3.0 Light, a PC program for backup and data base functions is included in the accessories.

MEP Instruments is sole distributor of leading names in sample analysis technology in Australia and New Zealand. The joint venture company formed by Metrohm Switzerland and Anton Paar Austria provides comprehensive service, support and training for the whole product line of Metrohm, Anton Paar, Autolab, Camag, Ismatec, Physica and Retsch in Australia and New Zealand. Contact MEP Instruments Pty Ltd (Australia), Unit 5, 112-118 Talavera Road, North Ryde, NSW 2113 Tel 02 9878 6900 Fax 02 9878 6999 [info@mep-instruments.com.au](mailto:info@mep-instruments.com.au); or in New Zealand: MEP Instruments Ltd (New Zealand) Brooke House, 24 Balfour Road, Parnell, Auckland Tel (+64) 09 366 1236 Fax (+64) 09 366 1235 [info@mep-instruments.co.nz](mailto:info@mep-instruments.co.nz) Visit MEP online: <http://www.mep.net.au>