

MEP Instruments Media Release

For immediate publication

Media Contact:

Sinikka D teindre, MEP Instruments Pty Ltd Sydney 02 9878 6900, sd@mep-instruments.com.au
Daniel Ledergerber, Creatvie Spiders Sydney, 02 9590 4777, ledi@creativespiders.com.au

Technical Contact Metrohm Titration:

Carlo Chiara, MEP Instruments Pty Ltd Victoria, 03 9810 4516, cc@mep-instruments.com.au

Photo: Metrohm KFT Titrino 795

Number of words: 238

Media release (incl. pictures in print resolution) can be downloaded from
<http://www.mep.net.au/press/release>

795 KFT Titrino for Karl Fischer titrations with real-time curve representation

The 795 KFT Titrino is a high-class compact KF Titrator with an LCD screen. Parameters can be viewed and altered "live" during the titration. All the functions are programmed via the separate keyboard, which can be removed for routine work. Operation is then carried out using the three keys on the front panel.

And additionally...

Apart from the well-known comfortable Titrino features such as simple programming of operator-specific methods, statistical evaluation, recalculations with altered parameters, the possibility of using sample changers, etc. the 795 KFT Titrino offers you the following additional possibilities:



- Method memory for over 100 operator methods.
- Sample data storage for automatic sample processing using sample changers with up to 156 sample data records.
- Alphanumerical input of method names, sample identifiers and result units or descriptions.
- Connection of a 776/765 Dosimat for automatic titer determination or automated back titrations.
- Connection for KF oven and 774 Oven Sample Processor.
- Barcode reader and PC keyboard can be connected (options).
- Printout conforming to GLP with result, titer, date and instrument serial number.

PC connection

The two built-in RS 232C interfaces are used for connecting a balance, a printer or a PC for additional data processing, result storage and data export to LIMS systems.

Metrodata VESUV' 3.0 Light, a PC program for backup and data base functions is included in the accessories.

MEP Instruments is sole distributor of leading names in sample analysis technology in Australia and New Zealand. The joint venture company formed by Metrohm Switzerland and Anton Paar Austria provides comprehensive service, support and training for the whole product line of Metrohm, Anton Paar, Autolab, Camag, Ismatec, Physica and Retsch in Australia and New Zealand. Contact MEP Instruments Pty Ltd (Australia), Unit 5, 112-118 Talavera Road, North Ryde, NSW 2113 Tel 02 9878 6900 Fax 02 9878 6999 info@mep-instruments.com.au; or in New Zealand: MEP Instruments Ltd (New Zealand) Brooke House, 24 Balfour Road, Parnell, Auckland Tel (+64) 09 366 1236 Fax (+64) 09 366 1235 info@mep-instruments.co.nz Visit MEP online: <http://www.mep.net.au>