

MEP Instruments Media Release

For immediate publication

Media Contact:

Sinikka D teindre, MEP Instruments Pty Ltd Sydney 02 9878 6900, sd@mep-instruments.com.au
Daniel Ledergerber, Creatvie Spiders Sydney, 02 9590 4777, ledi@creativespiders.com.au

Technical Contact Metrohm Titration:

Carlo Chiara, MEP Instruments Pty Ltd Victoria, 03 9810 4516, cc@mep-instruments.com.au

Photo: Metrohm Titrando 809_809

Number of words: 336


Media release (incl. pictures in print resolution) can be downloaded from
<http://www.mep.net.au/press/release>

808/809 Titrando Systems — experience the magic touch in titration

The basic unit consists of a Titrando with one measuring interface. This basic unit can be extended to form a super-titrator that controls 12 burets and boasts 6 galvanically separated measuring interfaces. Do you work already with Metrohm Exchange and Dosing Units equipped for automatic valve switching? No problem, you still have the possibility of using existing Exchange and Dosing Units models as well as 685 Dosimats and 700 Dosinos.

Touch control or mouse click?

No matter whether you operate the Titrando with the Touch Control and your fingers or via a PC and a mouse — both appearance and operation are practically identical. Methods and determinations are 100 % compatible. The operator guidance system can be regarded as exemplary.



Quality and compatibility

GLP, GMP, 21 CFR Part 11 are becoming increasingly important in day-to-day laboratory life. The Titrando System is entirely oriented toward quality management in the laboratory and offers the following possibilities: Control of rights of access using login and password, protection of electronic records against accidental or intentional alteration and complete traceability. If it has been programmed accordingly the instrument will remind you about any validation or service work that is due.

Intelligent dosing la carte

The Dosing and Exchange Units set new standards with regards to reliability. Intelligence in the form of an inconspicuous Data Chip makes it possible. This chip is present in every new Exchange and Dosing Unit. The Titrando automatically reads from it all the data (type of reagent, titer, last titer determination and much more) that it needs to carry out the titration properly.

Data management is knowledge management

The data that your Titrando System provides you with can be printed out and stored in a conventional manner. However, all advanced data management

functions are available. This includes storage on a PCMCIA card, which provides a link between the PC world and any other instrument that are equipped correspondingly. Do you want to know more about the Titrandos?

www.titrando.com

MEP Instruments is sole distributor of leading names in sample analysis technology in Australia and New Zealand. The joint venture company formed by Metrohm Switzerland and Anton Paar Austria provides comprehensive service, support and training for the whole product line of Metrohm, Anton Paar, Autolab, Camag, Ismatec, Physica and Retsch in Australia and New Zealand. Contact MEP Instruments Pty Ltd (Australia), Unit 5, 112-118 Talavera Road, North Ryde, NSW 2113 Tel 02 9878 6900 Fax 02 9878 6999 info@mep-instruments.com.au; or in New Zealand: MEP Instruments Ltd (New Zealand) Brooke House, 24 Balfour Road, Parnell, Auckland Tel (+64) 09 366 1236 Fax (+64) 09 366 1235 info@mep-instruments.co.nz Visit MEP online: <http://www.mep.net.au>