

MEP Instruments Media Release

For immediate publication

Media Contact:

Sinikka D teindre, MEP Instruments Pty Ltd Sydney 02 9878 6900, sd@mep-instruments.com.au
Daniel Ledergerber, Creatvie Spiders Sydney, 02 9590 4777, ledi@creativespiders.com.au

Technical Contact Metrohm Titration:

Carlo Chiara, MEP Instruments Pty Ltd Victoria, 03 9810 4516, cc@mep-instruments.com.au

Photo: Metrohm Titrando 831

Number of words: 167

Media release (incl. pictures in print resolution) can be downloaded from
<http://www.mep.net.au/press/release>

Water determination right down to trace levels

The microprocessor-controlled low cost 831 KF Coulometer continues Metrohm's established line of Karl Fischer instruments for the determination of low-level water contents in liquids, solids and gases. No titer determination required. Seven dialogue languages are standard, as is the backlit screen with its graphics capabilities. The screen not only provides a clear presentation of the parameters, but also shows graphically the course of the KF determination as g water versus time.

Convincing features:

- Cell with or without diaphragm — the choice is yours
- Internal method memory for 100 methods
- Compatible with ISO 9000 and GLP
- One press of a button for renewing the KF solution (option)
- Connection possibility for a barcode reader and a PC keyboard
- Possibility of setting limits for the sample weight and for the result
- The 768 KF Oven offers additional possibilities (option)



Metrodata VESUV' 3.0 Light, a PC program for backup and data base functions is included in the accessories.

MEP Instruments is sole distributor of leading names in sample analysis technology in Australia and New Zealand. The joint venture company formed by Metrohm Switzerland and Anton Paar Austria provides comprehensive service, support and training for the whole product line of

Metrohm, Anton Paar, Autolab, Camag, Ismatec, Physica and Retsch in Australia and New Zealand. Contact MEP Instruments Pty Ltd (Australia), Unit 5, 112-118 Talavera Road, North Ryde, NSW 2113 Tel 02 9878 6900 Fax 02 9878 6999 info@mep-instruments.com.au; or in New Zealand: MEP Instruments Ltd (New Zealand) Brooke House, 24 Balfour Road, Parnell, Auckland Tel (+64) 09 366 1236 Fax (+64) 09 366 1235 info@mep-instruments.co.nz Visit MEP online: <http://www.mep.net.au>